INFRABED & BIOBLANKET PRICES 26 MAY 21



InfraBed Posture~Pillow



\$440

InfraBed Posture~Pillow for managing neck traction to help correct forward head posture. Slide up and over to heat-treat. help restore neck curve and to open shoulders. Turn it over or around for deep therapy into your neck and shoulders. Sleep on your back or side.

InfraBed Mini Pad square

\$440

\$1360

\$330



InfraBed Mini-Pad is a 45 cm square and robustly suits any flat applications like seats, beds and under foot. Light, strong and easy to move and place. Relieves aches. Helps healing. Sits horizontally, or vertically, for the back, front or sides of body. "Use it with prostate." (PK)

InfraBed Medium Massage Size



InfraBed Medium Massage Bed size 150cm by 60cm. Big enough for full body infrared therapy. Small enough (13 kg) to move, or carry. (bag) Suits therapists and home users. The most popular InfarBed bed. Offer* = Free BioBlanket King Single creates Cocoon Effect.

Energy BioBlanket King Single



BioBlanket King Single is the perfect partner for InfraBed use. 200 cm by 120 cm. Designed to help recovery and rejuvenation. 5 InfraBed layers complete the Cocoon effect that protects both externally and reflects internally. Features 2 sides: Yin passive / Yang active. Zipup cotton cover.

BioBlanket Covers

from \$88



BioBlanket Covers are a part of BioBlanket's seven-layers of construction. Features both a Quilted (Yang) side, and a Sheeted (Yin) side, made from 100% quality cotton. Washable. **Zip-up** and so easy to change. King Single \$88. Double \$144.

InfraBed EndoBelt / FlexiPad



InfraBed EndoBelt / FlexiPad 52cm x 25cm. Super versatile. The FlexiPad has a blue stretch velcro belt to secure it around body / pelvis. Drape FlexiPad on, under, round your body, or on seat, or in bed. Set your time and temperature. Place against body. Use the cover or clothing.

InfraBed Chair Pod



InfraBed ChairPod is 45cm by 85cm. Made to flex through the middle, to suit any chair or bed angles. Suits flat on bed, and folded into chairs. Light enough to move, or carry easily. (bag) A favourite with young and old. "Especially good when the arthritis acts up." (SS)

InfraBed King Single Bed



InfraBed Big King Single Bed 164cm by 80cm. Utilises 20 kg of pure Amethyst crystals for ideal InfraBed hyperthermal therapy. Its 15 layers reduce emf and help you activate Core Body Temperature and balance bio-energies. "Deeper Rest, Recovery and Wellness" (PK)

Energy BioBlanket King Double



*On Special

\$1660



BioBlanket King Double is a big 240cm by 200cm. The ideal 'Cocoon' size for a completely reflective Bio-enclosure for external protection, and internal holistic bio-energetic feedback. Folds right over for one person, or under, or over, a double bed. Both sides do reflect fully.

ENAR Therapy System



ENAR Therapy System has 2 university research studies showing often 'dramatic and sustained results' for chronic (neck) pain, related disability, and general health. "ENAR works well with InfraBed and a BioBlanket in combination."

\$440

\$880

INFRABED THERAPY FEATURES 26 MAY 21




InfraBed Core Activation




Core Activation is the central therapeutic purpose* for the InfraBed range of hyperthermal therapy products. Used for working both locally or generally* to deeply penetrate for natural bio-activation with far infrared & subtle energies. * "feeling and moving better"

Holistic Hyper-Thermia



Holistic Hyperthermia is a most natural and popular holistic health practise. It includes Dry and Wet Saunas, Sun Bathing, Hot Springs, Hot Baths, Hot Lamps and Hyper Thermal Therapy InfraBeds. Controls are precise and easy. Add BioBlanket for a 'Cocoon'.

Amethyst Far Infrared Profile



InfraBed's Amethyst grade A crystals, when heated, emit a subtle and deeply penetrating Far Infrared emission profile, around 9.2 numens. This is the 'sweet spot' considered ideal in bio-energetic thermal therapy. Infrared is natural, light-based, and anti-inflammatory in effect.

IMulti Layers Construction



InfraBed 15 layer construction is designed to protect from external emissions, incl EMF. Also to produce and reflect both ideal subtle and Far Infrared **Bio-therapeutic Emissions.** Several of your Infrabed layers, including a negative ions layer, are also in Energy BioBlankets

Big 3 ET InfraBed Warranties




Enlightened Therapies takes its 'Big 3 Guarantee' : seriously. 1 is First Quality Manufacture 2 is Three Year Warranty and 3 is Complete Satisfaction so: ET guarantees you significant improvement of your issue or a full refund* (within 40 days*) Ph toll free ET 1300 305 107

InfraBed Hyper Thermal Therapy




InfraBed Hyper-Thermal Therapy gives you influence on your Core Body Temperature with InfraBed's Time and Temperature controllers. These you can set and forget from 30 - 70 degrees C Temp and Timer from 3 to 8 hours, depending on types of ACHE.

InfraBed ACHE Settings



InfraBed has a range of setting options for your therapy needs. A cute use high-temp low-time. C hronic use mid-temp & time. Holistic low-temp & long-time. *E* nergetic room-temperature ie body temperature for any length of time, for holistic bioenergy protection and calming.

InfraBed 'Sun' Power



InfraBeds emit a natural 'Sun-like' Far Infrared, lightbased deep thermal therapy. Viruses and cancers do not like heat. Good core body temperature activates and supports organs, circulation, metabolism and immune system functionality.

Build Quality Certifications



Your InfraBed range is Quality assured to CE standards and is certified by SAA approval for safety electrical standards in Australia and New Zealand. (240V, 288 Watt, Amps 1.22A) Guaranteed against any faults, in manufacture or materials with the ET full 3 Year Warranty.

ET : Our Purpose



Colleen & Paul Keetley began ET as therapeutic device manufacturers with a purpose 'Helping Every Body Heal.' Non-invasive natural healing works this way so ET's goal is To Empower People to Safely Treat themselves and others.* NB Consult your medical doctor

Empowerment